



CFD-00111-G
January 6, 2020



HL7 Conformance Statement

Synapse VNA Software Version 7.X

Copyright©2014–2020
FUJIFILM Medical Systems U.S.A., Inc.
All rights reserved.

REVISION HISTORY

| Date | Revision | Revisions |
|-----------------|----------|---|
| 6-January-2020 | G | <ul style="list-style-type: none"> • Synapse VNA version 7.0. • Updated copyright. |
| 16-May-2019 | F | <ul style="list-style-type: none"> • Added FHIR support. |
| 4-January-2019 | E | <ul style="list-style-type: none"> • Updated copyright year. |
| 10-May-2018 | D | <ul style="list-style-type: none"> • Added Patient Worklist Messaging. • Updated copyright year. |
| 9-March-2017 | C | <ul style="list-style-type: none"> • Updated copyright year. |
| 5-December-2016 | B | <ul style="list-style-type: none"> • Updated to reflect current operation. <ul style="list-style-type: none"> ○ Support for HL7 version 2.5 ○ Single and multi-byte character set ○ Ideographic and phonetic name support ○ File Upload Messaging ○ OBX File Upload ○ Non-DICOM Availability Notification |
| 23-July-2015 | A | <ul style="list-style-type: none"> • V6.0.4 General Availability release. • Rebranded with FUJIFILM logo, disclaimer and contact information. |
| 27-May-2015 | NA | <ul style="list-style-type: none"> • V6.0.3 General Availability release. |
| 19-April-2015 | NA | <ul style="list-style-type: none"> • V6.0.2 General Availability release. |
| 07-January-2015 | NA | <ul style="list-style-type: none"> • Updated copyright year. • Added additional copyrights information. |
| 12-August-2014 | NA | GA release for 6.0. |



Table of Contents

| | |
|---|-----------|
| Revision History | 2 |
| Introduction..... | 5 |
| Purpose | 5 |
| Acronyms and Abbreviations..... | 5 |
| Inbound Messaging Overview | 6 |
| Patient Demographics Messaging | 6 |
| Patient Merge Messaging | 6 |
| Patient Worklist Messaging | 6 |
| Patient Change ID Messaging | 7 |
| Orders Messaging | 7 |
| SR Results Messaging..... | 7 |
| File Upload Messaging | 7 |
| Inbound Segment Mappings..... | 8 |
| MSH Segment | 8 |
| PID Segment..... | 8 |
| MRG Segment | 9 |
| OBR Segment | 9 |
| ORC Segment | 10 |
| ZDS Segment | 10 |
| OBX Segment (when used for SR Report)..... | 10 |
| OBX Segment (when used for File Upload) | 10 |
| Outbound Messaging Overview..... | 11 |
| Study Availability Notification Messaging | 11 |
| EMR Viewer Study Availability Notification..... | 11 |
| Non-DICOM Availability Notification | 11 |
| Image Delete Messaging..... | 11 |
| Study Verification Messaging | 12 |
| Outbound Segment Mappings..... | 13 |
| MSH Segment | 13 |
| PID Segment..... | 14 |
| OBR Segment | 14 |
| OBX Segment (for Study Availability Notification) | 14 |
| OBX Segment (for EMR Viewer Study Availability Notification)..... | 15 |
| OBX Segment (for Image Delete Notification)..... | 15 |
| QRD Segment (for Study Verification) | 16 |

| | |
|---|-----------|
| Outbound Messaging Acknowledgement | 17 |
| MSA Segment..... | 17 |
| PID Segment..... | 17 |
| FHIR Capabilities | 18 |
| Capabilities Statement URL | 18 |
| Resources | 18 |
| DocumentReference..... | 18 |
| Binary..... | 18 |

Introduction

PURPOSE

This document describes the HL7 messages supported by Synapse VNA. It is intended for users involved with system integration and/or software design. We assume that the reader is familiar with the terminology and concepts that are used in HL7 and the IHE Technical Framework.

Synapse VNA supports HL7 versions 2.2, 2.3, 2.3.1, 2.4, 2.5. It also supports FHIR v3.0.1.

Although the use of this conformance statement in conjunction with the HL7 standard is intended to facilitate communication with Synapse VNA, it is not sufficient to guarantee, by itself, the inter-operation of the connection between a 3rd party system (such as a HIS/RIS or EMR) and Synapse VNA.

ACRONYMS AND ABBREVIATIONS

| Item | Description |
|------|---|
| ADT | Admission, Discharge and Transfer |
| EVN | Event Type segment |
| FHIR | Fast Health Interoperability Resources |
| HIS | Hospital Information System |
| HL7 | Health Level 7 |
| MRG | Merge Patient Information segment |
| MSH | Message Header segment |
| OBR | Observation Request segment |
| OBX | Observation/Result segment |
| PMS | Practice Management System |
| PID | Patient Identifier, also HL7 patient identifier segment |
| ORM | Order Request message |



Inbound Messaging Overview

Synapse VNA listens for messages generated by the ADT system of record, typically the EMR, HIS, departmental information system or a PMS. It processes the message, and issues an (original mode) acknowledgement.

PATIENT DEMOGRAPHICS MESSAGING

Synapse VNA handles the following demographic messages without any additional customization required. Any set of messages with a valid PID segment can be selected for patient demographics messaging. Processing of each is handled identically:

- ADT^A01 – Patient Admit Notification
- ADT^A02 – Patient Transferred
- ADT^A03 – Patient Discharged
- ADT^A04 – Patient Registered
- ADT^A05 – Patient Pre-Admission

These messages can be configured to work identically as above or to not change patient identity (patient ID or MPI).

- ADT^A08 – Update Patient Information
- ADT^A28 – Add Patient Information

Only the MSH and PID segments of these messages are processed, all other segments are ignored.

PATIENT MERGE MESSAGING

Synapse VNA processes the following patient merge messages:

- ADT^A18 – Merge Patient Information
- ADT^A34 – Merge Patients
- ADT^A40 – Merge Patients

Only the MSH, PID and MRG segments of this message are processed, all other segments are ignored.

PATIENT WORKLIST MESSAGING

The Patient Worklist Messaging function is the standard operation to create patient worklist entries for Connex Mobile. This is done using ADT^A14 messaging. A worklist entry is created using the Patient Account Number (PID-18) as the encounter number. If the Patient Account Number is not found, no worklist entry is created.

The following is the default ADT Worklist message Synapse VNA supports:

- ADT^A14 –Pending Admit – creates a Patient Worklist entry and a patient if one does not exist
- ADT^A03 –Discharge Patient –removes the Patient Worklist entry; does not remove the patient. It creates a patient if no patient exists.

Patient Worklist Messaging uses the MSH and PID segments: all other segments are ignored.

PATIENT CHANGE ID MESSAGING

- The MRG segment identifies the patient to update, and the PID segment contains the new MPI or Patient ID value. ADT^A46 –Change MPI ID
- ADT^A47 – Change Patient ID

Only the MSH, PID and MRG segments of this message are processed, all other segments are ignored.

ORDERS MESSAGING

Synapse VNA can listen for HL7 order messages to change the accession number of studies and optionally validate received DICOM studies against matching orders.

- ORM^O01 – General Order Message

Only the MSH, PID, ORC, OBR, and ZDS segments of these messages are processed. All other segments are ignored. The ORC, OBR segments can repeat as long as they are always in pairs.

SR RESULTS MESSAGING

A results message conveys the interpretation of an imaging study. Internally, Synapse VNA converts the interpretation (or diagnostic report) into a DICOM Structured Report (SR).

- ORU^R01 – Observation Results

Resulting messaging differs for all other HL7 messaging in that it creates a DICOM SR Image.

FILE UPLOAD MESSAGING

This is similar to an SR message, except the OBX segment includes the URL or http location of a file to upload. See the OBX segment (File Upload) from the *HL7 Messaging Specifications*.

Inbound Segment Mappings

MSH SEGMENT

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|-----------------------|-----------------|---|
| 1 | Field Separator | Yes | Used by Synapse VNA as field separator. |
| 2 | Encoding Characters | Yes | Used by Synapse VNA as encoding characters. |
| 3 | Sending Application | Optional | Used for MSH-5 in ACK. |
| 4 | Sending Facility | Optional | Used for MSH-6 in ACK, used to look up HL7SourceClass object. |
| 5 | Receiving Application | Optional | Used for MSH-3 in ACK |
| 6 | Receiving Facility | Optional | Used for MSH-4 in ACK |
| 9 | Message Type | Yes | Determines type of message. Anything not processed is simply ACK'ed. |
| 10 | Message Control ID | Yes | Used for MSH-10, MSA-2 in ACK. |
| 11 | Processing ID | Yes | Used for MSH-11 in ACK. |
| 12 | Version ID | Yes | 2.2, 2.3, 2.3.1, 2.4, or 2.5. |
| 18 | Character set | Optional | Specifies single byte or multi-byte character sets found by scanning the HL7 Message for ISO 2022 escape sequences. See tables in the <i>HL7 Messaging Specifications</i> for supported character sets. |

PID SEGMENT

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|--------------------------|-----------------|--|
| 2 | Patient ID (External ID) | Optional | If present, is used as the patient's master patient index (MPI). |
| 3 | Patient ID (Internal ID) | Yes | Supplies the patient ID. Multiple patient ID's are allowed (~). |

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|--------------------------|-----------------|--|
| 5 | Patient Name | Yes | <p>First five subfields of name are stored in LastName, FirstName, MiddleName, Suffix, and Prefix.</p> <p>If patient already exists, patient name is always assumed to be complete. Null values in last four subfields will overwrite existing strings.</p> <p>For HL7 v2.5 only, we support multi-byte character encodings, which enables ideographic and phonetic writing systems.</p> |
| 7 | Date/Time of Birth | Optional | Stored to database. |
| 8 | Sex | Optional | Stored to database. |
| 10 | Race | Optional | Stored to database. |
| 18 | Patient Encounter Number | Optional | Stored to database. |

MRG SEGMENT

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|-----------------------------|-----------------|---|
| 1 | Prior Patient ID - Internal | Yes | Supplies the merge from patient ID. Multiple patient IDs are allowed (~). |
| 4 | Prior Patient ID - External | Optional | If present, is used as the source patient's master patient index (MPI). |

OBR SEGMENT

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|-----------------------|-----------------|---|
| 2 | Placer Order Number | Optional | Stored in database. |
| 3 | Filler Order Number + | Conditional | <p>Stored in database as Accession Number</p> <p>Note: While not required per the HL7 standard, this is required by Synapse VNA for proper operation.</p> |
| 4 | Universal Service ID | Yes | Stored in database. |
| 6 | Requested Date/Time | Optional | Stored in database. |

ORC SEGMENT

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|----------------|-----------------|--|
| 1 | Order Control | Yes | NW (Create order), CA (Cancel order), or XO (change study and order accession number). |

ZDS SEGMENT

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|----------------|-----------------|--|
| 1 | Study UID | Yes | Used to lookup study for XO operation. |

OBX SEGMENT (WHEN USED FOR SR REPORT)

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|---------------------------|-----------------|---|
| 1 | Set ID - OBX | Optional | Incrementing value starting at 1, one segment for SR results line. |
| 5 | Observation Value | Optional | Set to SR text line. Also can use tilde (~) to represent different lines. |
| 11 | Observation Result Status | Optional | Set to A if study is complete, C if study is updated. |

OBX SEGMENT (WHEN USED FOR FILE UPLOAD)

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|-------------------|-----------------|--|
| 1 | Set ID - OBX | Optional | Incrementing value starting at 1, one segment for SR results line. |
| 2 | Value Type | Optional | For upload messaging, we only process OBX segments with the type RP. |
| 5 | Observation Value | Optional | Filename to upload. Name is either in HTTP notation (string starts with http:// or https://) or UNC notation. |

Outbound Messaging Overview

Synapse VNA generates internal events for changes to its archives. These events can be used to trigger a custom outbound HL7 message:

- Study notification
- EMR notification
- Non-DICOM EMR notification
- Delete notification
- Study verification

STUDY AVAILABILITY NOTIFICATION MESSAGING

This message is sent to notify an external system that a study is available:

- ORU^R01 – Observation Results

The following segments are supplied: MSH, PID and OBR and one or more OBX.

EMR VIEWER STUDY AVAILABILITY NOTIFICATION

Synapse VNA can send a URL link to a study's DICOM images or non-DICOM files when they are available for viewing:

- ORU^R01 – Observation Results

The following segments are supplied: MSH, PID and OBR and one or more OBX.

NON-DICOM AVAILABILITY NOTIFICATION

This is based on the EMR Viewer Study Availability Notification. The URL takes you to the study in the non-DICOM mode.

IMAGE DELETE MESSAGING

This message is sent whenever a study, series, or image(s) are deleted:

- ORU^R01 – Observation Results

The following segments are supplied: MSH, PID and OBR and one or more OBX.

STUDY VERIFICATION MESSAGING

Another way that Synapse VNA interacts with an external Enterprise system is to verify patient and study information. If the verification fails, Synapse VNA marks the study with a QA issue so an administrator can later resolve the discrepancy between the data Synapse VNA received from the DICOM modality and the Enterprise's data:

- OSQ^Q06 – Order Status Query

The following segments are supplied: MSH and QRD.

The response is expected to be an OSR^Q06 with the MSH, MSA, QRD, and PID segments present (which will be parsed).

Outbound Segment Mappings

MSH SEGMENT

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|---------------------------------|-----------------|---|
| 1 | Field Separator | Yes | Used by Synapse VNA as field separator. |
| 2 | Encoding Characters | Yes | Used by Synapse VNA as encoding characters. |
| 3 | Sending Application | Optional | From SendingApplication attribute in the HL7Client object. |
| 4 | Sending Facility | Optional | From SendingApplication attribute in the HL7Client object. |
| 5 | Receiving Application | Optional | From Receiving Application attribute in the HL7Client object. |
| 6 | Receiving Facility | Optional | From Receiving Application attribute in the HL7Client object. |
| 9 | Message Type | Yes | Set to ORU^R01. |
| 10 | Message Control ID | Yes | Set to incrementing number. |
| 11 | Processing ID | Yes | Set to P. |
| 12 | Version ID | Yes | 2.2, 2.3, 2.3.1, 2.4, or 2.5 |
| 15 | Accept Acknowledgment Type | Optional | Set to AL. |
| 16 | Application Acknowledgment Type | Optional | Set to NE. |

PID SEGMENT

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|--------------------------|-----------------|---|
| 2 | Patient ID (External ID) | Optional | From the database. |
| 3 | Patient ID (Internal ID) | Yes | From the database. |
| 5 | Patient Name | Yes | From the database. First five subfields of name come from in LastName, FirstName, MiddleName, Suffix, and Prefix. For HL7 v2.5 only, we support multi-byte character encodings, which enables ideographic and phonetic writing systems. |
| 7 | Date/Time of Birth | Optional | From the database. |
| 8 | Sex | Optional | From the database. |
| 10 | Race | Optional | From the database. |

OBR SEGMENT

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|-------------------------|-----------------|--------------------|
| 1 | Set ID--OBR | Optional | Set to 1. |
| 2 | Placer Order Number | Optional | From the database. |
| 4 | Universal Service ID | Yes | From the database. |
| 7 | Observation Date/Time # | Conditional | From the database. |
| 24 | Diagnostic Serv Sect ID | Optional | From the database. |

OBX SEGMENT (FOR STUDY AVAILABILITY NOTIFICATION)

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|------------------------|-----------------|-----------------------------------|
| 1 | Set ID - OBX | Optional | Incrementing value starting at 1. |
| 2 | Value Type | Conditional | Set to RP. |
| 3 | Observation Identifier | Yes | From the database. |

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|------------------------------|-----------------|---|
| 4 | Observation Sub-ID | Conditional | Incrementing value for study, starting at one. |
| 5 | Observation Value | Optional | From the database. |
| 6 | Units | Optional | Set to number of images. |
| 7 | References Range | Optional | Set to series description. |
| 11 | Observation Result Status | Yes | Set to A if study is complete, C if study is updated, D if notification is for deleted images. |
| 13 | User Defined Access Checks | Optional | Set to ONLINE if study stored in online storage group. Set to NEARLINE if study stored only in nearline storage group. |
| 14 | Date/Time of the Observation | Optional | Set to series created time. |
| 15 | Producer's ID | Optional | Set to series station name. |
| 17 | Observation Method | Optional | Set to series modality. |

OBX SEGMENT (FOR EMR VIEWER STUDY AVAILABILITY NOTIFICATION)

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|-------------------|-----------------|--|
| 1 | Set ID - OBX | Optional | Set to 1. |
| 2 | Value Type | Conditional | Set to RP if configured to send EMR URL, otherwise set to TX . |
| 5 | Observation Value | Optional | Set to EMR Viewer URL if configured, otherwise 'ONLINE' . |

OBX SEGMENT (FOR IMAGE DELETE NOTIFICATION)

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|-------------------|-----------------|-------------------------------------|
| 1 | Set ID - OBX | Yes | Incrementing value starting at one. |
| 2 | Value Type | Conditional | Set to RP . |
| 5 | Observation Value | Conditional | From the database. |

QRD SEGMENT (FOR STUDY VERIFICATION)

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|---------------------------|-----------------|---|
| 2 | Query Format Code | Optional | Set to R. |
| 3 | Query Priority | Optional | Set to I. |
| 4 | Query ID | Optional | Set to 1. |
| 8 | Who Subject Filter | Yes | Set to patient to be validated. |
| 7 | Date/Time of Message | Optional | Set to current date/time. |
| 10 | What Department Data Code | Optional | Set to study to be validated, using accession number. |

Outbound Messaging Acknowledgement

The external information systems should respond with the following message segments:

MSA SEGMENT

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|----------------------|-----------------|---|
| 1 | Acknowledgement Code | Yes | If set to anything but AA , verification failed. |
| 2 | Message Control ID | Optional | If failure, text recorded with QA issue. |
| 3 | Text Message | Optional | If failure, text recorded with QA issue. |

PID SEGMENT

| SEQ | HL7 Field Name | Required by HL7 | Synapse VNA Usage |
|-----|--------------------|-----------------|---|
| 2 | Patient ID | Optional | Stored to database. |
| 5 | Patient Name | Optional | First five subfields of name are stored in LastName, FirstName, MiddleName, Suffix, and Prefix. |
| 7 | Date/Time of Birth | Optional | Stored to database. |
| 8 | Sec | Optional | Stored to database. |
| 10 | Race | Optional | Stored to database. |

FHIR Capabilities

Fast Health Interoperability Resources (FHIR) is a standard for healthcare data exchange. It is a global standard published by HL7 International based on RESTful web services. Synapse VNA acts as a FHIR server to enable non-DICOM image storage and retrieval. This section provides a summary of those capabilities. Full details are available from the VNA Capabilities Statement URL below.

FHIR sources are linked to an organization within the VNA. For authentication, FHIR uses the same security model as the non-DICOM RESTful web services use. FHIR supports the HTTP/HTTPS protocol, with HTTPS being the default.

CAPABILITIES STATEMENT URL

<VNA Base URL>/fhir/metadata

RESOURCES

DOCUMENT REFERENCE

Interactions supported:

- Create
- Read
- Update
- Delete

BINARY

Interactions supported:

- Read