

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
FUJIFILM Medical Systems U.S.A., Inc.	Flash IIP Console	Flex UI v3.1	1/20/2006
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow	Acquisition Modality	N/A	
Patient Information Reconciliation			
Basic Security (ATNA)	Secure Node		
<b>Internet address for vendor's IHE information:</b> <a href="http://www.fujimed.com/sub/product_bulletins.html">http://www.fujimed.com/sub/product_bulletins.html</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>			
<b>DICOM</b>	<a href="http://www.fujimed.com/sub/product_bulletins.html">http://www.fujimed.com/sub/product_bulletins.html</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	