

# SYNAPSE

CV 6.0 and above

## IHE Integration Statement

\* SYNAPSE is a registered trademark of FUJIFILM MEDICAL SYSTEMS, U.S.A., Inc. in United States, Australia and Canada.

\* SYNAPSE is a registered trademark of FUJIFILM Corporation in China, Japan, Austria, Belgium, Denmark, Finland, France, Germany, Greece, Ireland, Italy, Luxembourg, The Netherlands, Portugal, Spain, Sweden, United Kingdom, Latvia, Slovakia, Slovenia, Cyprus, Lithuania, Hungary, Czech Republic, Malta, Poland, Estonia, Bulgaria, & Romania.

---

**Manufacturer:**

FUJIFILM CORPORATION  
26-30, Nishiazabu 2-Chome Minato-Ku,  
Tokyo 106-8620 JAPAN

---

**Legal Manufacturer's Address:**

FUJIFILM Medical Systems U.S.A., Inc.  
81 Hartwell Avenue, Suite 300  
Lexington, Massachusetts 02421 U.S.A.

---

**Manufacturing Site:**

FUJIFILM Medical Systems U.S.A., Inc.  
3020 Carrington Mill Blvd, Suite 500  
Morrisville, North Carolina 27560 U.S.A.

---

**European Authorized Representative:**

FUJIFILM Europe GmbH  
Heesenstrasse 31  
40549 Düsseldorf, Germany

---

**Notified Body:**

Intertek Semko AB  
Torshamnsgatan 43, Box 1103  
SE-164 22 Kista, Sweden

---

**Global Support:**

<http://www.fujifilm.com/worldwide/>

---

# C E0413

Copyright © 1997-2019 FUJIFILM MEDICAL SYSTEMS, U.S.A., Inc. All rights reserved.

Unauthorized reproduction or distribution of this program, or any portion of it, is strictly prohibited.

Other products and company names mentioned herein may be trademarks under the respective owners.

---

**Proprietary Information:** For Distribution within FUJIFILM Medical Systems U.S.A., Inc. and its Authorized Affiliates, Distributors, Partners, and Customers Only.  
Duplication or Distribution of this document without agreement from FUJIFILM Medical Systems U.S.A., Inc. is strictly prohibited.



Synapse IHE Integration Statement

DOC-0042770-A

Page 3 of 3

IHE Integration Statement

Vendor	Product Name	Version	Date
FUJIFILM	Synapse CardioVascular	6 and above	06/17/2019
This product is intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Image Manager / Archive Image Display Performed Procedure Step Manager Department System Scheduler/Order Filler	None	
Patient Information Reconciliation	Image Manager / Archive Performed Procedure Step Manager Department System Scheduler/Order Filler	None	
Cardiac Cath Workflow	Image Manager / Archive Image Display Performed Procedure Step Manager Department System Scheduler/Order Filler	None	
Echocardiography Workflow	Image Manager / Archive Image Display Performed Procedure Step Manager Department System Scheduler/Order Filler	None	
Retrieve ECG for Display	Information Source Display	None	
<b>Internet address for vendor's IHE information:</b> <a href="http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/">http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/</a>			
<b>Links to Standards Conformance Statements of the Implementation</b>			
DICOM	<a href="http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/">http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/</a>		
HL7	<a href="http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/">http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/</a>		
Links to general information on IHE			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">www.ihe-j.org</a>	